

Compact Attract

GO ANYWHERE, DO ANYTHING



The küschall Compact attract meets the expectations of users who are looking for a reliable, cross-foldable active wheelchair, offering all the adjustability required for secure seating, maneuverability and easy driving.



Yes, you can.®

küschall®

Compact Attract



The Compact attract makes daily tasks really simple. Our proven quality is reflected in this wheelchair, and it is the perfect companion for your activities. Ingenious ideas make it compact. Its precision makes it robust.

Features



Transport size by a small width.



Vertical adapter plate for easy adjustment of the rear wheel position



T-armrest with soft pad and clothes guards.



Push brake in open position.

Tech Specs



Seat width
40, 45, 50cm



Seat depth
40, 45, 50cm



Seat height
front: 430 - 520 mm
rear: 400 - 490 mm



Total width
seat width + 18cm



Total width folded
30cm



Total length
70°: approx. 985 mm
90°: approx. 895 mm
w/out legrest:
690 mm



Max. user weight
130 kg



Transport weight
8.7 kg



Total weight unloaded
11.3 kg



Turning radius
92cm



Max. safe slope
7°



Crash Test Approved
Fulfills the requirements of ISO 7176-19



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au

Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com

Thailand: Phone +66 0 2108 8555 or visit www.invacare.com



Yes, you can.®