

# We do Beds



**ALEGIO NG**

**CS7**

**CS7 KING SINGLE**

**ETUDE PLUS**

<b>Bed Length:</b>	2110mm	2140mm	2140mm	2190mm
<b>Bed Width</b>	900mm	900mm	1090mm	1020mm
<b>Adjustable Height:</b>	280 - 800mm	170 - 760mm	170 - 760mm	330 - 730mm or 400 - 800mm
<b>Max. User Weight:</b>	135kg	226kg	226kg	165kg
<b>Mattress Length:</b>	1970mm	1970mm	1970mm	1970mm
<b>Mattress Width:</b>	880mm	880mm	1070mm	880mm
<b>Bed Weight</b>	78kg	86kg	92kg	77kg



**ETUDE PLUS LOW**

**OCTAVE**

**1600**

**SCANBETA NG**

<b>Bed Length:</b>	2190mm	2221mm	2060mm	1680 or 1880mm
<b>Bed Width</b>	1020mm	1190mm	920mm	780 or 880mm
<b>Adjustable Height:</b>	200 - 600mm or 270 - 670mm	400 - 800mm	200 - 670mm	300 - 900mm
<b>Max. User Weight:</b>	165kg	350kg	165kg	70kg
<b>Mattress Length:</b>	1970mm	2000mm	1970mm	1600 or 1800mm
<b>Mattress Width:</b>	880mm	1160mm	880mm	700 or 800mm
<b>Bed Weight</b>	77kg	170kg	170kg	91kg or 99.5kg

For more information on the Invacare Beds range call **0800 468 222** or visit [www.invacare.co.nz](http://www.invacare.co.nz)



**Yes, you can.®**

**Making Life's Experiences Possible**