

Ocean Ergo Accessories



ACCESSORIES | ADJUSTMENTS | COMPATIBILITY



Yes, you can.®

AQUATEC®

Seating and backrests



Standard soft seat
With key shaped cut out.



Soft seat insert
Can be used with or without the soft seat when commode function is not required.



Small soft seat
With small, oval opening for smaller and young adults.



Small soft seat insert
Can be used with the small soft seat when the commode function is not required.



Ergonomic soft seat
Extra soft, textured surface for pelvic stability.



Variable soft seat
Hygiene cut out can be positioned in all 4 directions.



Universal soft seat
Comfortable, closed seat surface.



XL back
Increases the width between the armrests by 80mm. Can be retrofitted to standard Ocean Ergo models.



Solid plastic backrest
Offers a stable upright position and is easy to clean.



Wider XL armrests
Increases the distance between the armrests by 40 mm on each side.



Soft backrest cushion
For additional stability and comfort. Will also benefit smaller users and young adults.



Comfort backrest
For additional comfort and easy cleaning.

Safety



Safety bar
For additional safety and reassurance.



Safety belts Ocean
For safe positioning. Individually adjustable, high load capacity and easy to assemble. Machine washable.



Armrest lock
Locks armrest into the lower or higher position.



Front and rear anti-tippers
For additional safety.

Positioning



Lateral supports

For safe positioning. Width and height adjustable – also for small users and young adults.



Headrest holder

Safe and comfortable with various adjustment positions (pad not included).



Ergonomic headrest/neckrest

Offers exceptional management of the head position. Neckrest shown.



Calf support with foot support

Padded, both height and angle adjustable with additional foot support.

Toileting



Calf and amputee support

Padded, both height and angle adjustable. For use with the standard plastic seat, standard soft seat and small soft seat.



Splash guard **NEW DESIGN**

Clips on easily for use with the standard plastic seat, standard soft seat and small soft seat.



Toilet pan **NEW DESIGN**

- Ergonomic handle allows carer to remove pan from holder easily
- Scale inside the pan
- Lid holder when lid is not required.

Adjustments

All Ocean Ergo models can be adjusted to suit a range of user sizes.

ARMRESTS

Ergonomic, swing away armrests can be adjusted to different heights.

FOOTRESTS

Height adjustable, swing away and detachable footrests. No tools required.

SEAT HEIGHT

Seat to ground height can be adjusted between 475 mm - 600 mm. No tools required.

FOOTPLATES

Angle adjustable, flip-up and textured footplates with heel loops. No tools required.



Yes, you can.®

Compatibility

Most accessories are cross compatible with all Ocean Ergo models.	Product Code	Ocean Ergo	Ocean 24" Ergo	Ocean Ergo XL	Ocean 24" Ergo XL	Ocean VIP Ergo	Ocean Dual VIP Ergo
Standard soft seat	AP1603142	✓	✓	✓	✓	✓	✓
Soft seat insert	AP1603143	✓	✓	✓	✓	✓	✓
Small soft seat	AP1603145	✓	✓	✓	✓	✓	✓
Small soft seat insert	AP1603152	✓	✓	✓	✓	✓	✓
Splash guard	AP1603155	✓	✓	✓	✓	✓	✓
Universal soft seat	AP1603154	✓	✓	✓	✓	✓	✓
Ergonomic special soft seat	1531259	✓	✓	✓	✓	✓	✓
Variable soft seat	1544469	✓	✓	✓	✓	✓	✓
Safety bar	10167-10	✓	✓	X	X	✓	✓
Safety bar XL	1470078	X	X	✓	✓	X	X
Chest belt	1470081	✓	✓	✓	✓	✓	✓
Pelvic belt	1470082	✓	✓	✓	✓	✓	✓
Armrest locking device	1470077	✓	✓	✓	✓	✓	✓
XL Back	AP1604015	✓	✓	-	-	X	X
Fixed backrest	1575299	✓	✓	X	X	✓	✓
Soft backrest cushion	1558769	✓	✓	✓	✓	✓	✓
Wider XL armrests	1470086	✓	✓	X	X	✓	✓
24" wheels (set)	16379	✓	-	X	-	X	X
Calf and amputee support right	AP1643573	✓	✓	✓	✓	✓	✓
Calf and amputee support left	AP1643574	✓	✓	✓	✓	✓	✓
Calf support with foot support right	AP1643543	✓	✓	✓	✓	✓	✓
Calf support with foot support left	AP1643572	✓	✓	✓	✓	✓	✓
Lateral	1535077	✓	✓	✓	✓	✓	✓
Ergonomic headrest holder	1531613	X	X	X	X	✓	✓
Ergonomic headrest	1531465	X	X	X	X	✓	✓
Ergonomic neckrest	1531608	X	X	X	X	✓	✓
Toilet pan	10230	✓	✓	✓	✓	✓	✓
Toilet pan holder	AP1602078	✓	✓	✓	✓	✓	✓
Front and back anti-tipper Ergo VIP models	AP1603992	X	X	X	X	✓	✓
Front anti-tipper Ocean Ergo models	AP1601986	✓	✓	✓	✓	X	X
Comfort soft backrest VIP Ergo models	AP1603606	X	X	X	X	✓	✓
Comfort soft backrest Ocean Ergo models	AP1604136	✓	✓	X	X	X	X
Comfort soft backrest Ocean Ergo XL models	AP1604016	X	X	✓	✓	X	X

✓ Compatible X Not Compatible - Included



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz
Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au
Asia: Phone +61 2 8839 5330, email anzexportsales@invacare.com or visit www.invacare.com
Thailand: Phone +66 0 2821 5515 or visit www.invacare.co.th



Yes, you can.®