




# We do Mobility Scooters

## Colour Options:



	COLIBRI	LEO	PEGASUS METRO	PEGASUS PRO	COMET PRO (NZ ONLY)	COMET ALPINE+	COMET ULTRA
<b>Backrest Width:</b>	465mm	470mm	510mm	510mm	510mm	510mm	660mm
<b>Seat Depth:</b>	400mm	410mm	470mm	470mm	470mm	460mm	535mm
<b>Floor to Seat Height:</b>	370 - 420mm	460 - 560mm	440 - 510mm	440 - 510mm	440 - 510mm	440 - 510mm	440 - 510mm
<b>Backrest Height:</b>	385mm	475mm	490mm	490mm	490mm	630mm	515mm
<b>Unit Driving Width:</b>	505mm	590mm	660mm	650mm	665 - 685mm	665 - 685mm	665 - 685mm
<b>Total Length:</b>	1010mm	1220mm	1270mm	1320mm	1450mm	1450mm	1450mm
<b>Empty Weight:</b>	49kg (18 Ah)	59kg (exc. batteries) 83kg (inc. batteries)	110kg	136kg	136kg	143kg	148kg
<b>Max. User Weight:</b>	136kg	136kg	136kg	160kg	160kg	160kg	220kg
<b>Battery Capacity:</b>	2 x 12 V/18Ah	2 x 12V/36Ah	2 x 12 V/50Ah	2 x 12 V/75Ah	2 x 12 V/75Ah	2 x 12 V/75Ah	2 x 12 V/75Ah
<b>Turning Diameter:</b>	2200mm	2620mm	2600mm	2750mm	2750mm	2750mm	2750mm
<b>Max. Climable Obstacle Height:</b>	45mm	60mm	80mm	100mm	100mm	100mm	100mm
<b>Ground Clearance:</b>	50mm	70mm	78mm	100mm	100mm	100mm	100mm
<b>Max Safe Slope:</b>	6°	10°	6°	10°	10°	10°	10°
<b>Max. Driving Distance*:</b>	Up to 16km	38km	54km	52km	55km	58km	49km
<b>Max. Driving Speed:</b>	8 km/h	8 km/h	10 km/h	10 km/h	15 km/h	10 km/h	10 km/h
<b>Motor Capacity:</b>	200 W	240 W	240W / 600W	550 W / 1300 W	550 W / 1600 W	650 W / 1800 W	650 W / 1800 W

\* Maximum driving distance is dependent on driving habits, terrain and battery age.



For further information and full technical specifications please see your dealer or visit [www.invacare.co.nz](http://www.invacare.co.nz) or [www.invacare.com.au](http://www.invacare.com.au)



**Making Life's Experiences Possible**