

Küschall® K-SERIES



| UNLIMIT YOUR WORLD



UNLIMIT YOUR WORLD

We at Küschall believe anything is possible when you have the power to achieve it. The Küschall K-Series offers modularity and performance in a single product. Unlimit your Strength, Unlimit your Freedom, Unlimit your world.

KÜSCHALL K-SERIES

The K-series is an icon of design, a trendsetter in the monotube frame construction. The ingenious modular seating system will give you the perfect fit and its dynamic set up will offer you a superb driving performance, to keep you on the move.

HYDROFOR- MING

The hydroformed front frame will bring the design of your K-Series to another level. The new technology brings more stiffness and rigidity, for an unrivaled driving performance, adding striking lines and a clear design language to a timeless icon.



Küschall®

K-SERIES



Yes, you can.®

Küschall®

K-SERIES

HYDROFORMED KNEE AREA

- Ergonomic grip, ideal for repositioning and transfers, frame angle available in both 90° and 75°.

LIGHT ACTIVE BRAKE

- The new aluminum brake, with a minimalist design, brings lightness to unprecedented levels. Precision even in the smallest details.

TAPERED FRAME

- Tapered frame allows the wheelchair to contour to your body for a sleek aesthetic.

FOOTBOW

- Completely redesigned with a flat surface for greater comfort while driving, and perfect positioning.





NEW BACKBRACE

- Ergonomic handle, a crucial aspect in the chair lifting and unloading phases. Lower brace positioning, ideal for Matrix backrests.

BACKREST

- Folding backrest with six possible angle adjustments for different positioning needs.

SEAT MODULE

- At the heart of the K-series. The seat module offers five different positions of the centre of gravity, making the wheelchair more active.

CASTOR HOUSING

- New aluminum profile welded to the front frame tube. Consistent design throughout the whole wheelchair.

FORK

- New clean design and integrated bearings. Collecting dirt will significantly be reduced and maintenance will be easier.

KEY FEATURES

ICONIC

With its iconic minimalist open frame and lightweight design, the K-Series gives you the ability to be able to enjoy everyday activities; whether it's getting to and from work, or pursuing your exciting hobbies, the K-Series gives you the freedom to be you.



MODULAR

An invisible centre of gravity adjustment beneath the seat means you can move the seat module back and forth into five different positions. It can be easily modified at a later date with no extra components required.

DYNAMIC

The frame has a constant perfect geometry to the ground, so there is no need to adjust anything at all on the frame.

Wheels always remain in the same position on the frame. Very rigid construction, leads to great driving performance.




HYDROFORMED FRONT FRAME

Hydroforming technology has added even more personality to the K-Series. An Iconic look that today becomes even more recognizable thanks to the straight edge design of the knee area. With the frame further strengthened, your driving experience will be more thrilling than ever.

CUSTOMISATION

1 CHOOSE THE FRAME COLOR



Pure White



Anthracite Grey



Space Black



Azure Blue



Carmine Red

2 CHOOSE THE DECAL STYLE



3 CHOOSE THE DECAL COLOR



Black



Red



Light Grey



Gold



Red



Blue


















Dark Grey



Black

TECH SPECS

<p>SEAT WIDTH</p> <p>320 - 500 mm</p> <p>In increments of 20 mm</p> 	<p>SEAT DEPTH</p> <p>350 - 525 mm</p> <p>In increments of 25 mm</p> 	<p>LOWER LEG LENGTH</p> <p>290 - 480 mm</p> <p>in increments of 10 mm</p> 	<p>BACKREST HEIGHT</p> <p>270 - 480 mm</p> <p>In increments of 15 mm</p> 	<p>BACKREST ANGLE</p> <p>74° / 78° / 82° / 86° / 90° / 94°</p> 
<p>SEAT HEIGHT</p> <p>front: 450 - 520 mm</p> <p>rear: 390 - 490 mm</p> <p>in increments of 10 mm</p> 	<p>FRAME ANGLE</p> <p>75° / 90°</p> <p>V-frame or tapered</p> 	<p>TOTAL WIDTH</p> <p>Seat width + 170 mm</p> 	<p>TOTAL LENGTH</p> <p>75°: approx. 830 - 1055 mm</p> <p>90°: approx. 760 - 985 mm</p> 	<p>CENTER OF GRAVITY</p> <p>80 mm / 95 mm / 110 mm / 125 mm / 140 mm</p> 
<p>CAMBER</p> <p>0° / 3° / 6°</p> 	<p>TOTAL WEIGHT UNLOADED</p> <p>approx. 7.8 kg</p> 	<p>TRANSPORT WEIGHT</p> <p>approx. 5.1 kg</p> 	<p>MAX. USER WEIGHT</p> <p>130 kg</p> 	<p>CRASH TEST</p> <p>Approved - Fulfills the requirements of ISO 7176-19</p> 

All measurements and weights indicated are based on seat width 400 mm and seat depth 400 mm in the lightest and most compact configuration.

Weight and measurements can change with the different configurations. For further information please check the user manual on our website.



Yes, you can.®

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication.
Invacare reserves the right to alter product specifications without prior consent. 05/2024



Yes, you can.[®]