

Action 1R

A FEATURE RICH, ENTRY LEVEL CHAIR



An entry level wheelchair that doesn't compromise on functionality.



Yes, you can.®

Making Life's Experiences Possible



Action 1R

The Action 1R is an entry level wheelchair that doesn't compromise on functionality. Built using decades of in-house design experience, the Action 1R is available in transit and self-propelling versions in three sizes. The removable, swing away armrests are easy to remove and the armpad can be used in two positions.

Features



With years of experience and using tried-and-tested technology, the Action 1R is a reliable partner, the one you can count on.



The optional Eco elevating legrests provide extra comfort for users.



Optional cane holder can be mounted to securely hold a cane or crutches.



Dynamic anti-tippers provide extra safety on slopes or when backrest is reclined.

Tech Specs



Seat height
49.5cm



Seat depth
42.5cm



Total width
SW + 20cm
(self-propel)
SW + 16cm (transit)



Total length
108cm



Seat width
40.5, 45.5, 50.5cm



Max. user weight
125kg



Backrest height
43cm

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 05/2024



Yes, you can.®