

Solara 3G

TILT-IN-SPACE WHEELCHAIR



The Invacare® Solara® 3G Wheelchair retains this durability, and once again revolutionizes tilt-in-space wheelchair design with incredible simplicity, portability and ease of use and adjustment.



Yes, you can.®

Making Life's
Experiences Possible



Solara 3G

Features

- 136kg standard, 180kg weight capacity (heavy-duty package)
- Center-of-gravity remains constant throughout tilt range
- Positive lock mechanism eliminates “sponginess” common to locking hydraulic cylinders
- Tilt angle indicator
- Standard flared front provides more space for user’s feet and reduces interference when tilting in the anterior direction
- One axle plate accommodates all wheel sizes
- -5° to 50° of total tilt available
- Adjustable seat depth, seat width and back cane angle.

Options



The **DeGage Dynamic Rocker Back** accommodates more aggressive rockers, like children with uncontrolled extension of the upper and lower extremities. The Rocker Back will allow an individual to move in the direction they are pushing and as a result will help to reduce pain or discomfort. For more information please see www.degage.us.

Tech Specs



Seat width
12 - 24"



Seat depth
12 - 22"



Seat height
12.5 - 19"



Tilt-In-Space range
55°



Total weight unloaded
40kg



Turning Radius
25"



Max. user weight
136kg
181kg (HD)



Max. safe slope
9°



Crash Test Approved
Fulfills the requirements of ISO 7176-19

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 05/2024



Yes, you can.®