

Action[®] Ampla[™]

PROVIDING INDEPENDENCE FOR AN ACTIVE LIFESTYLE



**TAKING BARIATRIC WHEELCHAIRS TO
THE NEXT LEVEL**

INNOVATIVE | AGILE | AWARD WINNING





Taking bariatric wheelchairs to the next level

The Invacare **Action® Ampla™** has been designed and developed specifically for Plus-Size users, taking into consideration body shapes, varying activity levels and weight limit requirements. But it doesn't stop there! The sleek and minimalist wheelchair, with its modern appearance, provides independence and support for a comfortable and active lifestyle.

Award winning

Already an award-winning chair, the agile and stylish Invacare **Action® Ampla™** has won a Red Dot Award for Product Design. This award is not only the largest product competition in the world, but also one of the most prestigious design awards. The **Action® Ampla™** has also won the Italian Platinum A'Design Award and the Gold German Design Award.





Yes, you can.®

5

KEY FEATURES

1

Inbuilt transportation solution

The transportation bearings on the **Action® Ampla™** have been cleverly designed to take the weight of the wheelchair, allowing it to easily go through doorways or narrow spaces:

1. Fold down the backrest
2. Lift the legrest section of the wheelchair
3. Tilt wheelchair onto the transportation bearings

Whilst in this position, the wheelchair can be moved easily.

1.



2.



3.



5

Multiple active rear wheel positions

- ▶ Very easy to propel
- ▶ Users can go further, using less energy



2

Ergonomic backrest

- ▶ Maximum access to handrims
- ▶ No obstructions
- ▶ More freedom of movement
- ▶ Superb driving performance and manoeuvrability



3

Height adjustable transfer levers

- ▶ Support for entering and leaving the chair
- ▶ Can be turned inwards or outwards
- ▶ Can be lowered to allow table access



4

Personalised accent

- ▶ Three stylish colour choices





Yes, you can.®

Innovative Backrest Design



The design of the wheelchair allows for a perfect fit for both Plus-Size male users, who tend to have more of an "apple" body-shape and Plus-Size female users, who tend to have a "pear" body-shape. This is thanks to the innovative backrest on the **Action® Ampla™** that has been designed to accommodate the anatomy of both body-shapes whilst offering support and comfort to the user.



Main features













- ▶ Inbuilt transportation solution
- ▶ Ergonomic backrest
- ▶ Height adjustable transfer lever
- ▶ Personalised accent
- ▶ Multiple active rear wheel positions

Additional features


- ▶ An elastic, breathable privacy cover hides the gluteal shelf
- ▶ Push to lock brakes are nicely integrated into the side guard for accessibility

Tech Specs


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

 <p>Seat width 50.5cm, 55.5 cm, 60.5cm, 65.5cm, 70.5cm, 75.5cm</p>	 <p>Seat depth 46-51cm</p>	 <p>Seat height 38 - 45.5cm in steps of 2.5cm</p>	 <p>Backrest height 45 - 52.5cm</p>	 <p>Total width seat width + 20cm (self-propel) seat width + 14.5cm (transit)</p>	 <p>Total length 108cm</p>
 <p>Max. user weight 300kg</p>	 <p>Transport weight From 19.5kg</p>	 <p>Total weight From 27kg</p>	 <p>Turning radius 142 - 166cm</p>	 <p>Max. safe slope 7°</p>	 <p>Crash Test Approved Fulfills the requirements of ISO 7176-19</p>

Frame Colours



Accent Colours



Please note: colours may vary slightly from those displayed here.



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication.
Invacare reserves the right to alter product specifications without prior consent. **05/2024**



Yes, you can.®