

e-fix

ENJOY A NEW DEGREE OF MOBILITY



Turn your manual wheelchair into a compact electric wheelchair with the e-fix add-on drive. The electric add-on drive is the simple way to increase your freedom of movement and can be individually programmed to your personal needs. The e-fix also allows you to cover greater distances and handle slopes effortlessly.



Yes, you can.®





e-fix

It is so easy to control your wheelchair with e-fix you can get used to it almost at once. Buttons that are clearly arranged and easy to reach on the control unit ensure the easiest possible handling. You always have all the important information clearly in view on the easily legible display. The backlit display can be read easily, even in the dark. The small control unit ensures more comfort when driving. The speed can be varied continuously through a controller that is easy to grip. The lower the speed selected, the more agile the driving behaviour. You steer and the e-fix thinks for you!

Features



Powerful Drive Wheels

Easy to fix onto the wheelchair with quick-lock release fasteners and just as easy to remove again in a moment to save space in transport.



Ergonomic Control Unit

With large and reliable battery status display, individually adjustable driving characteristics and ergonomically arranged controls.



Small, Light Battery Pack

Suspended under the seat to save space and allows a range of up to 16 km*.

* Depending on the type of wheelchair, weight of the occupant, type of ground, wheelbase and the centre of gravity of the wheelchair

Tech Specs

Battery pack	36 V
Top Speed	6 km/h
Climbing capability	up to 20 % restricted to wheelchair model
Wheel dimensions	22, 24"
max. user weight	120, 160 kg
max. permitted total weight	170, 210 kg (463 lb)
Weight of components:	
Battery pack	2 kg (4.4 lb)
Drive wheel	8.0 kg (17.6 lb)
Joystick control unit	0.6 kg (1.3 lb)

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 04/2024



Yes, you can.®