

e-pilot

ADD-ON SCOOTER FOR MANUAL
WHEELCHAIRS



The e-pilot turns your wheelchair into a sporty vehicle in a matter of seconds. Just a few simple steps to attach the e-pilot to your wheelchair and you have a dynamic combination for everyday life. The e-pilot offers an environmentally friendly mobility alternative with maximum driving pleasure!



Yes, you can.®





e-pilot

Take a quick trip into the city, visit friends or just a quick run with the dog: the e-pilot makes it possible. Even longer distances are not an issue for you due to the large lithium-ion battery pack. When you arrive at the destination, you simply undock the e-pilot and use your manual wheelchair as usual.

The e-pilot offers you a new kind of mobility. The large front wheel allows you to expand your options: A dirt road or cobblestones are no longer obstacles. Even over a curb, not a problem. The sturdy chassis, integrated battery pack and powerful drive wheel combine to make a sporty design with sophisticated ergonomics perfectly.

Features



Battery pack

Auto contacting, quick-release battery pack, range up to 50 km (depending on the ground surface, the weight of the user, temperature, wheelchair settings and terrain), weight: only 2.9 kg, battery can be charged in the installed or removed state, optionally available spare battery pack for even more range.



Automatic battery-charger

Electronic charging and switch-off mechanism, charging time of an empty battery approx. 4 hours.



Display

Easy to remove and attach by twist-lock, TFT full colour display with indicator for speed, battery capacity, driving mode and tour data, language selection, traction assist and system information via menu.



Colour options

Individual appearance through a wide range of design colours. Frame colour is black. Colour panels on the side frame and fork are white as standard or choose among a variety of colours. Colour choice for colour panels:
 - 4 glossy signal colours
 - 6 Matt metallic colours

Tech Specs

Max. User Weight	Max. 100kg
Max. Total Weight	Max. 135kg
Climbing Capability	Up to 10%
Braking Capability	Up to 20%
Battery	36 V
Top Speed	20 km/h
Range	Up to 50 km
Total Weight of e-Pilot	18.4kg

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 04/2024



Yes, you can.®