

Alegio NG

OPTIMISED COMFORT



The Alegio NG fits easily into any home care or institutional environment. The lowest position means that shorter or more unstable users can easily and safely get up, or out of, the bed. The maximum height, at 80 cm, ensures healthy working conditions for care staff.



Yes, you can.®

Making Life's
Experiences Possible

Alegio NG

- **Large under bed clearance** - Good for OBT's and Lifter access.
- **Optimised mattress sections** - Fits 99% of adult western population.
- **Compliant with IEC 60601-2-52** - Safety first for patient and carer peace of mind.
- **Tool-less assembly** - Easy to install by hand, with no loose parts and assemble and a quick release bracket between the base and the mattress platform.
- **Slat deck design** - Improved infection control and reduces wear and tear on mattresses.
- **Adjustable height range of 280 - 800mm** - The lowest position allows shorter patients to easily and safely get up, or out of the bed.



Features



Removable head and foot boards.



Two electrical function with manually operated leg section.



Tool-less assembly



Easy to store and to transport when placed on the transport kit.

Accessories








Folding metal side rail
Code: 1554518-0154



Lifting pole
Code: 1542093-0154

Tech Specs

				
Length	Width	Adjustable Height	Max. User weight	Bed weight
2110mm	900mm	280 to 800mm	135kgs,	78kg
1970mm	880mm			

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 04/2024



Yes, you can.®