

Dacapo Heavy Duty

At Risk

FIRM COMFORTABLE AND CUSHIONED

Designed to match customer's requirements within Invacare's homecare and hospital bed ranges. With a 350kg maximum user weight, the Dacapo® Heavy Duty has been developed to offer optimal cushioning and pressure relief without compromising on hygiene and user friendliness.



Unique Features

Benefits

Heavy duty	Tailored to meet the needs of the larger user and offers a firm, comfortable cushioned support
Heavy user cover	100% fluid-proof, coated polyester with polyurethane, Breathable, Two way stretch capacity, Washable at 95 degrees
Two high resilience foams	This model can be reversed to adjust to the degree of firmness and support required by the user
Fire retardancy	Conforming to European Standard EN 597 - 1 and European Standard EN 597 - 2

Features



Heavy Duty
maximum 350 kg:
Top (50 mm) visco-elastic foam 50 kg/m³:
Base (50 mm) high resilient foam 38 kg/m³:
Base (50 mm) high resilient foam 50 kg/m³



Easy Handling.




Discretely zippered on three sides with cloth covering.



Designed to fit the Octave bariatric bed.


Tech Specs

PRODUCT CODE: 1451259



Length

2050mm



Width

1160mm



Depth

150 mm



User weight

350kg

WARRANTY

Warranty
2yr cover
2yr foam

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 04/2024



Yes, you can.®