

MX2 Back

PEAK PERFORMANCE AND LOOKS GREAT



The Invacare® Matrx® MX2 Back is the latest in back support technology. Designed for active users who demand the best, it combines aesthetic appeal with superior support and function. Lightweight carbon fibre is specially engineered to provide unmatched stability for peak performance.



Yes, you can.®

matrx®
SEATING SERIES

MX2 Back



Cover:

- Breathable Spacer Mesh fabric, 9" and higher backs include leatherette side panels
- Cover and foam overlap back shell for privacy, protection and comfort
- Includes removable lumbar positioning pad
- Optional privacy flap available

Carbon Fiber Back Shell:

- Exceptional durability and aesthetic appeal
- Very light weight
- 3" Contour to provide stability to hips/trunk without restricting upper body mobility








Hardware:

- Choice of Fixed Mount or Quick Release Hardware. Hardware does not come standard with the MX2 back, this will need to be selected.

Accessories



Tech Specs

MX2 Back	 Height 6, 9, 12, 16"	 Width 14 - 18"	 Depth adjustment 2" Fixed HW 4" Quick-Release HW	 Angle adjustment 30°	 Product weight from 680g (1.5lb)	 Weight Capacity 114kg (250lb)	 Limited Warranty 2 years
-----------------	---	---	--	---	---	--	---

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 07/2024



Yes, you can.®