

PB Deep Back

EXCEPTIONAL SUPPORT AND COMFORT



The Invacare® Matrix® PB Deep Back is contoured to provide increased midline support and positioning. The extensive range of sizes provides exceptional comfort and support for a wide range of individuals with moderate to high lateral trunk positioning needs.



PB Deep Back



Foam:

- Molded foam wraps around shell for maximum comfort and protection
- Includes removable positioning sacral/lumbar support pad

Outer cover:

- Moisture resistant and breathable
- Optional privacy flap is available

Back shell:

- Contoured shell provides exceptional midline support and positioning
- Durable, lightweight ABS back shell includes 7" contour depth and built-in plate for headrest attachment

Hardware:

- "Easy Set" hardware installed and adjusted in minutes - adjustments can be made while user is seated in the wheelchair

Accessories



Tech Specs

PB Deep	 Height 12", 16", 20"	 Width 14", 16", 17", 18", 20"	 Depth adjustment 2" forwards, 2" backwards	 Angle adjustment 40° range**	 Product weight From 1.3kg (2.9lb)	 Weight Capacity 136kg (300lb)	 Limited Warranty 2 years
PB Deep Heavy Duty	16", 20"	20, 22, 24, 26, 28, 30"	2" forwards, 2" backwards	40° range**	from 2.7kg (6lb)	226kg (500lb)	2 years

** Limitation of depth/angle adjustment if used on chairs narrower than model width

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 05/2024



Yes, you can.®