

# Elite Back

DESIGNED WITH ACTIVE LIFESTYLES IN MIND



The Invacare® MatrX® Elite Back is designed to provide precise orientation within the wheelchair for optimal postural support. The attractive, lightweight design appeals to individuals with active lifestyles.



Yes, you can.®

**matrx**®  
SEATING SERIES

# Elite Back



## Outer cover:

- Oversized cover/foam for extra comfort and protection
- Moisture resistant and breathable
- Includes removable lumbar positioning pad
- Optional privacy flap available

## Back shell:

- Stylish, lightweight and breathable aluminium back shell with 3" contour depth.

## Hardware:

- "Easy Set" hardware installed and adjusted in minutes - adjustments can be made while user is seated in the wheelchair

## Accessories



## Tech Specs

<b>Elite</b>	 <b>Height</b> 10", 12", 14"*, 16", 18", 20"	 <b>Width</b> 12" - 20"	 <b>Depth adjustment</b> 2" forwards, 2" backwards	 <b>Angle adjustment</b> 40° range***	 <b>Product weight</b> from 800g (1.7lb)	 <b>Weight Capacity</b> 136kg (300lb)	 <b>Limited Warranty</b> 2 years
<b>HD Elite Heavy Duty</b>	<b>Height</b> 16" - 20"	<b>Width</b> 21" - 24"	<b>Depth adjustment</b> 2" forwards, 2" backwards	<b>Angle adjustment</b> 40° range***	<b>Product weight</b> From 2.7kg (6lb)	<b>Weight Capacity</b> 227kg (500lb)	<b>Limited Warranty</b> 2 years

\* 14" height only available up to 18" wide

\*\*\* Limitation of depth/angle adjustment if used on chairs narrower than model width

FOR FURTHER INFORMATION:

**New Zealand:** Freephone 0800 468 222, email [sales@invacare.co.nz](mailto:sales@invacare.co.nz) or visit [www.invacare.co.nz](http://www.invacare.co.nz)

**Australia:** Freephone 1800 460 460, email [orders@invacare.com.au](mailto:orders@invacare.com.au) or visit [www.invacare.com.au](http://www.invacare.com.au)



All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 09/2024

**Yes, you can.®**