

Elite Deep Back

DESIGNED WITH ACTIVE LIFESTYLES IN MIND



The Invacare® Matrix® Elite Deep Back is designed with integrated lateral supports to provide optimal trunk support while allowing clearance for the shoulders and hips. The attractive, lightweight design appeals to individuals with active lifestyles.



matrix[®]
SEATING SERIES

Elite Deep Back



Outer cover:

- Oversized cover/foam for extra comfort and protection
- Moisture resistant and breathable
- Includes removable lumbar positioning pad
- Optional privacy flap available

Back shell:

- Stylish, lightweight and breathable aluminium back shell with 6" contour depth

Hardware:

- "Easy Set" hardware installed and adjusted in minutes - adjustments can be made while user is seated in the wheelchair

Accessories



Tech Specs

Elite Deep Back	 Height 12", 16", 18", 20"	 Width 14" - 20"	 Depth adjustment 2" forwards, 2" backwards	 Angle adjustment 40° range***	 Product weight from 1.3kg (3lb)	 Weight Capacity 136kg (300lb)	 Limited Warranty 2 years
Elite Deep Back Heavy Duty	Height 16, 20"	Width 21 - 24"	Depth adjustment 2" forwards, 2" backwards	Angle adjustment 40° range***	Product weight From 2.7kg (6lb)	Weight Capacity 227kg (500lb)	Limited Warranty 2 years

* 14" height only available up to 18" wide

*** Limitation of depth/angle adjustment if used on chairs narrower than model width

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 05/2024



Yes, you can.®