

# Ocean Ergo

AN ADVANCED RANGE OF MODULAR SHOWER CHAIRS



## PREMIUM PERSONAL CARE SOLUTIONS

ERGONOMIC | VERSATILE  
EASY TO MANOEUVRE



Yes, you can.®

AQUATEC®

# Even greater versatility and ergonomic enhancements

Our latest Aquatec Ocean Ergo range has been specifically designed to make showering and toileting even easier.

From a standard attendant propelled shower chair, to self-propel or a tilt-in-space, Ocean Ergo offers a solution for almost everyone.

In collaboration with leading therapists, our design goals were focused on the ergonomics, lifestyle and wellbeing of users.

Every model in the Ocean Ergo range is optimised with new designs to achieve higher levels of safety, comfort and support.

Whether being used independently or with the help of a carer, each one of our highly modular chairs offers a dignified solution.

## Tilt-in-space models provide additional supportive seating solutions

The Ocean VIP Ergo and Ocean Dual VIP Ergo have a  $-5^{\circ}$  to  $40^{\circ}$  seat tilt angle with an integrated weight shift mechanism that ensures optimal stability on a small footprint.

The Ocean Dual VIP Ergo is multi-functional benefiting from both tilt and recline. This allows the seat and backrest angle to open up between  $5^{\circ}$  and  $35^{\circ}$  in a simple smooth adjustment thanks to the lightweight frame.



# Introducing our latest **Ocean** Ergo range

---

Discover the full range below including optional accessories for specific requirements.



## **Ocean Ergo**

An entry level attendant propelled model



## **Ocean Ergo XL**

An entry level attendant propelled model with a heavier maximum user weight



## **Ocean 24" Ergo**

An entry level self-propelling model



## **Ocean 24" Ergo XL**

An entry level self-propelling model with a heavier maximum user weight



## **Ocean VIP Ergo**

A tilt-in-space model for additional seating and positioning support



## **Ocean Dual VIP Ergo**

A tilt-in-space model with recline for additional seating and positioning support

## NEW features

### Improved armrest system

- > Newly designed armrests provide **increased support and comfort for the user**. The new and optional armrest pads offer individuals a softer surface area and resting position. Featuring an easy to operate optional integrated safety lock system.



### Reductions in seat height

- > Additional **seat to ground height adjustment** has been introduced to promote safe and independent transfers. The frame can be lowered further accommodating people of a smaller stature, available on the Ocean Ergo XL.



### Tilt stop function

- > An easy to apply optional function which **prevents the frame tilting anteriorly past 0°**. This ensures the frame remains horizontal to the ground and is an aid for those who are hoisted in/out of the chair.



### Ergonomic headrest holder

- > With multiple axis adjustment the ergonomic headrest holder offers **superior postural support for the head/neck region**. Adjustments are quick and easy by simply rotating the integrated hand-screws to achieve the correct position for the individual. Ideal for assessment purposes and can be set up with no tools required.





## Innovative ergonomic seat

- > **Insight-led design** encourages a more upright seating position for enhanced independence and comfort.
- > Designed around the **imprints** people leave behind when seated.
- > A **5° seat tilt** lifts the knees into a natural squat position for toileting.
- > Transfer supports at front to allow the user to position themselves comfortably.



## Highly configurable

- > Available with an extensive range of positioning, toileting and seating accessories that can be used to support a wide range of conditions and postural needs.



## Extremely stable

- > **Stainless steel frame** for enhanced stability, additional safety and reassurance.



## Easy to handle

- > A small, compact footprint and an **optimised wheel base** makes it easy to manoeuvre in small spaces.
- > Easy disassembly for transportation and storage.



Yes, you can.®

# Tech Specs

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.



Seat width ▼



Seat depth ▼



Total width ▼



Total depth ▼



Width between armrests ▼

	Seat width ▼	Seat depth ▼	Total width ▼	Total depth ▼	Width between armrests ▼
<b>Ocean Ergo</b>	480 mm	480 mm	565 mm	905 mm	475 mm
<b>Ocean Ergo XL</b>	480 mm	480 mm	645 mm	905 mm	555 mm
<b>Ocean 24" Ergo</b>	480 mm	480 mm	670 mm	1050 mm	475 mm
<b>Ocean 24" Ergo XL</b>	480 mm	480 mm	670 mm	1050 mm	555 mm
<b>Ocean VIP Ergo</b>	480 mm	480 mm	565 mm	1000 mm	475 mm
<b>Ocean Dual VIP Ergo</b>	480 mm	480 mm	565 mm	1000 mm	475 mm



Tilt angle seat ▼



Tilt angle back ▼



Weight ▼



Load capacity ▼

	Tilt angle seat ▼	Tilt angle back ▼	Weight ▼	Load capacity ▼
<b>Ocean Ergo</b>	-	-	17 kg	150 kg
<b>Ocean Ergo XL</b>	-	-	18 kg	180 kg
<b>Ocean 24" Ergo</b>	-	-	23 kg	150 kg
<b>Ocean 24" Ergo XL</b>	-	-	24 kg	180 kg
<b>Ocean VIP Ergo</b>	-5° - 40°	-	23.5 kg	150 kg
<b>Ocean Dual VIP Ergo</b>	-5° - 40°	5° - 35°	28.3 kg	150 kg



Frame colour



Yes, you can.®

FOR FURTHER INFORMATION:

**New Zealand:** Freephone 0800 468 222, email [sales@invacare.co.nz](mailto:sales@invacare.co.nz) or visit [www.invacare.co.nz](http://www.invacare.co.nz)

**Australia:** Freephone 1800 460 460, email [orders@invacare.com.au](mailto:orders@invacare.com.au) or visit [www.invacare.com.au](http://www.invacare.com.au)

All rights reserved. All information quoted is believed to be correct at time of publication.  
Invacare reserves the right to alter product specifications without prior consent. **04/2024**



**Yes, you can.®**