

# Comet Ultra

BUILT TO SUPPORT WITH PERFORMANCE  
AND RELIABILITY



Ideal for those with larger body shapes, the Comet® Ultra offers a secure, supportive and powerful driving experience without compromising on style or functionality. The Comet Ultra will get you to your destination quickly and efficiently, in complete comfort and control.



Yes, you can.®

Making Life's  
Experiences Possible

# Comet Ultra



The safety and stability of the Comet Ultra is never a concern; even for users weighing up to 220kg. Developed with the highest quality components and safety features, the Comet Ultra makes a trustworthy companion no matter where you go. Furthermore, its flexibility and adaptability means you can customise your scooter to suit your unique style and individual preference.

The special contoured Ultra Seat has been developed in conjunction with therapists to ensure an optimum seating position and increased comfort. In addition, a pneumatically adjustable lumbar support allows for further comfort when seated if and when required. The reinforced seat base and adjustable armrests offer excellent support and allows the use of the armrest as additional support when getting on and off the scooter.

## Features



LCD dashboard for clear, easy to read information.



Improved suspension to offer superior comfort to larger individuals.



















Reinforced seat column offers additional stability and safety when seated.




Reinforced height and width adjustable armrest.

## Tech Specs

 <b>Backrest width</b> 660 mm	 <b>Seat depth</b> 535 mm	 <b>Floor to seat height</b> Min. 440 mm Max. 510 mm	 <b>Backrest height</b> 515 mm	 <b>Driving unit width</b> 665 - 685 mm	 <b>Total length</b> 1450 mm	 <b>Empty weight</b> 148 kg	 <b>Max user weight</b> 220 kg
 <b>Battery capacity</b> 2 x 12 V/75 Ah	 <b>Turning diameter</b> 2750 mm	 <b>Max climbable obstacle height</b> 100 mm	 <b>Ground clearance</b> 100 mm	 <b>Max safe slope</b> 10°	 <b>Max driving distance</b> 49 km	 <b>Max driving speed</b> 10 km/h	 <b>Motor capacity</b> 650 W / 1800 W

## Colour Options



\* Warranty Terms and Conditions apply. Contact Invacare for full details.

FOR FURTHER INFORMATION:

**New Zealand:** Freephone 0800 468 222, email [sales@invacare.co.nz](mailto:sales@invacare.co.nz) or visit [www.invacare.co.nz](http://www.invacare.co.nz)

**Australia:** Freephone 1800 460 460, email [orders@invacare.com.au](mailto:orders@invacare.com.au) or visit [www.invacare.com.au](http://www.invacare.com.au)

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 03/2025



**Yes, you can.®**