



TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Post Code _____ Country _____

Contact Name _____ Comments _____

	Code	Description	Price
<input type="radio"/>	CS1567978	Action ³ NG 40.5x40.5cm, Transit <i>Equipped with 12" transit wheels inc. fixed axles</i>	
<input type="radio"/>	CS1567975	Action ³ NG 40.5x40.5cm, Self-propelled <i>Equipped with 24" self propelling wheels inc. quick-release axles</i>	
<input type="radio"/>	CS1552784	Action ³ NG 40.5x45cm, Transit <i>Equipped with 12" transit wheels inc. fixed axles</i>	
<input type="radio"/>	CS1552711	Action ³ NG 40.5x45cm, Self-propelled <i>Equipped with 24" self propelling wheels inc. quick-release axles</i>	
<input type="radio"/>	CS1567979	Action ³ NG 45x40.5cm, Transit <i>Equipped with 12" transit wheels inc. fixed axles</i>	
<input type="radio"/>	CS1567976	Action ³ NG 45x40.5cm, Self-propelled <i>Equipped with 24" self propelling wheels inc. quick-release axles</i>	
<input type="radio"/>	CS1552843	Action ³ NG 45x45cm, Transit <i>Equipped with 12" transit wheels inc. fixed axles</i>	
<input type="radio"/>	CS1552769	Action ³ NG 45x45cm, Self-propelled <i>Equipped with 24" self propelling wheels inc. quick-release axles</i>	
<input type="radio"/>	CS1567980	Action ³ NG 50.5x40.5cm, Transit <i>Equipped with 12" transit wheels inc. fixed axles</i>	
<input type="radio"/>	CS1567977	Action ³ NG 50.5x40.5cm, Self-propelled <i>Equipped with 24" self propelling wheels inc. quick-release axles</i>	
<input type="radio"/>	CS1552844	Action ³ NG 50.5x45cm, Transit <i>Equipped with 12" transit wheels inc. fixed axles</i>	
<input type="radio"/>	CS1552770	Action ³ NG 50.5x45cm, Self-propelled <i>Equipped with 24" self propelling wheels inc. quick-release axles</i>	
FRAME INCLUSIONS			
	CEA 0980	Anthracite Black colour	N/C
	CEA 1010	8" x 1 1/4" castors, solid	N/C
	CEA 1270	Push to lock brakes	N/C
SEAT INCLUSIONS			
	CEA 0340	Nylon upholstery with velcro strips	N/C
	CEA 0270	48.5cm seat height	N/C
BACK INCLUSIONS			
	CEA 0410	Fixed 7° back	N/C
	CEA 0475	Lumbar support (+1 strap) upholstery, padded shaped cover	N/C
LEGREST INCLUSIONS			
	CEA 0510	80° straight swing in/out hanger RH	N/C
	CEA 0560	80° straight swing in/out hanger LH	N/C
	CEA 0600	Flip-up footplate RH	N/C
	CEA 0620	Flip-up footplate LH	N/C
ARMREST INCLUSIONS			
	CEA 0690	Swing-away detachable, height-adjustable full-length armrests RH <i>Adjustable 17-29.5cm</i>	N/C
	CEA 0790	Swing-away detachable, height-adjustable full-length armrests LH <i>Adjustable 17-29.5cm</i>	N/C
SAFETY INCLUSIONS			
	CEA 1360	Buckle belt, mounted	N/C
	CEA 1370	Rear anti-tippers, pair	N/C
	CEA 1385	Rear tie-down hooks	N/C
OPTIONAL EXTRAS			
ELEVATING LEGRESTS			
<input type="radio"/>	Angle adj legrest with angle/depth adj footplate RH	<input type="radio"/> 40cm (AP1543228R) <input type="radio"/> 45cm (AP1543230R) <input type="radio"/> 50cm (AP1543231R)	
<input type="radio"/>	Angle adj legrest with angle/depth adj footplate LH	<input type="radio"/> 40cm (AP1543228L) <input type="radio"/> 45cm (AP1543230L) <input type="radio"/> 50cm (AP1543231L)	
<input type="radio"/>	Amputee legrest RH (1442990)		
<input type="radio"/>	Amputee legrest LH (1442991)		
TRAY			
<input type="radio"/>	Tray table	<input type="radio"/> 40cm (1450922) <input type="radio"/> 45cm (1450924) <input type="radio"/> 50cm (1450926)	
HEADREST			
<input type="radio"/>	Neckrest inc. backbrace	<input type="radio"/> 40cm (AP1554576) <input type="radio"/> 45cm (AP1554578) <input type="radio"/> 50cm (AP1554580)	
CANE HOLDER			
<input type="radio"/>	Cane holder (0572-02056)		



Print Form

Email Form

Clear Form

ISO, MDA, TUV Logos

The product includes all standard features

Maximum user weight: 125kg

Crossbar Folding Mechanism, Seat Angle Possibility: Delivered Always 3°

Revised July 2014

