

Rovi X3

THE COMPLETE SOLUTION



The ROVI X3 Power Base is the perfect complement to Invacare's full line of innovative power positioning products! The innovative and exciting ROVI X3 boasts the narrowest wheelbase in the industry.



Rovi X3




- Industry's most modular system, the Ultra Low Maxx Seating System, featuring up to 55° of CG Tilt, and up to 12" of nearly silent Power Seat Elevation
- Available with the full range of Motion Concepts' accessories and Invacare® Matrix® Seating products
- Up to 136kg weight capacity
- Seat-to-Floor Heights as low as 17" with 12" of Power Seat Elevation
- 11" - 22" of seat width and 12" - 24" of seat depth adjustment*
- All powered seating functions can be operated either through or separate from the driver control.

Tech Specs

| | | | | | |
|--|-------------------------------------|------------------------------------|---------------------------------|---------------------------------------|--------------------------------------|
| Seat width 12 - 22" | Seat depth 12 - 24" | Seat height 17 - 20" | Base length 36.5" | Base width 23.25" | Max. user weight 136kg |
| Battery capacity Group 30 (60Ah) | Max. safe slope 8 deg | Motor type 4 pole | Top Speed 9km/h | Total weight from 130kg | Turning radius 20" |

Base Shroud Colours



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au

Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com

Thailand: Phone +66 0 2108 8555 or visit www.invacare.com



Yes, you can.®