

The Invacare Bed Range:



	ALEGIO	CS8 FLOORLINE SINGLE	CS7	ETUDE PLUS
Bed Length:	2110mm	2110mm	2140mm	2190mm
Bed Width	900mm	920mm	900mm	1020mm
Adjustable Height:	300 - 850mm	150 - 680mm	170 - 760mm	330 - 730mm or 400 - 800mm
Max. User Weight:	135kg	250kg	226kg	165kg
Mattress Length:	2000mm	1980mm	2030mm	2000mm
Mattress Width:	900mm	860mm	900mm	900mm



	ETUDE PLUS LOW	OCTAVE	1600	SCANBETA NG
Bed Length:	2190mm	2221mm	2060mm	1680 or 1880mm
Bed Width	1020mm	1190mm	920mm	780 or 880mm
Adjustable Height:	200 - 600mm or 270 - 670mm	400 - 800mm	200 - 670mm	300 - 900mm
Max. User Weight:	165kg	350kg	200kg	70kg
Mattress Length:	2000mm	2000mm	1890mm	1600 or 1800mm
Mattress Width:	900mm	1160mm	880mm	700 or 800mm

For more information on the Invacare Beds range call **1800 468 222** or visit www.invacare.com.au



Yes, you can.®

Making Life's Experiences Possible