



# Rea Azalea Minor

SOLID AND STABLE

The Rea Azalea Minor is designed for people that require a seat that has a smaller depth and width. Adapted from the Rea Azalea, the Azalea Minor boasts all the advantages of a reliable, tilt-able wheelchair with its special weight-shifting mechanism.







# Rea Azalea Minor

Rea Azalea Minor is designed for people that need a smaller-sized seat, such as petite adults or teenagers.

One of the key features of Azalea Minor is the 'Tilt in Space' weight-shifting mechanism, which ensures a stable tilt and makes manoeuvring easier. As with all the Azalea models, the width and depth of the Azalea Minor can be adjusted to suit individual needs.

## Features



Table tray is transparent and the rail offers gripping opportunities for posture solutions.



The chassis is 30 mm shorter, the seat depth is 50 mm shorter and the backrest 100 mm lower than on standard models.

## Tech Specs



**Seat width**  
34 - 39cm



**Seat depth**  
38 - 45cm



**Seat height**  
40cm



**Total width**  
seat width + 22cm



**Total length**  
97cm



**Legrest angle**  
Adjustable



**Tilt Range**  
25°



**Max. user weight**  
75kg



**Transport weight**  
14.5kg



**Total weight unloaded**  
32kg



**Crash Test Approved**  
Fulfills the requirements of ISO 7176-19



FOR FURTHER INFORMATION:

**New Zealand:** Freephone 0800 468 222, email [sales@invacare.co.nz](mailto:sales@invacare.co.nz) or visit [www.invacare.co.nz](http://www.invacare.co.nz)

**Australia:** Freephone 1800 460 460, email [sales@invacare.com.au](mailto:sales@invacare.com.au) or visit [www.invacare.com.au](http://www.invacare.com.au)

**Asia:** Phone +61 2 8839 5330, email [exportsales@invacare.com.au](mailto:exportsales@invacare.com.au) or visit [www.invacare.com](http://www.invacare.com)

**Thailand:** Phone +66 0 2108 8555 or visit [www.invacare.com](http://www.invacare.com)



**Yes, you can.®**