

ROVI | X3 Maxx



The complete solution

The all new **ROVI X3 Power Base** is the perfect complement to Motion Concepts' full line of innovative power positioning products! Available exclusively through Motion Concepts, the innovative and exciting **ROVI X3** boasts the narrowest wheelbase in the industry and all the function and durability you have come to expect from Motion Concepts.



Yes, you can.


Ultra Low Maxx System

- Industry's most modular system, the Ultra Low Maxx Seating System, featuring up to 55° of CG Tilt, 170° of Recline with Extended Shear Reduction (ESR) and up to 12" of nearly silent Power Seat Elevation
- Available with the full range of Motion Concepts' accessories and Invacare® Matrx® Seating products, mounted via our all new Ultra Rail and Ultra Track slotted interfaces
- Up to 181kg weight capacity
- Seat-to-Floor Heights as low as 17" with 12" of Power Seat Elevation
- 16" - 24" of seat width and 15" - 22" of seat depth adjustment
- All powered seating functions can be operated either through or separate from the driver control.



Also available - Tilt/Recline/Elevate, Elevate 12" and 6".

Technical Data




**Rehab seat and 70° swingaway footrest

Base Shroud Colour*

Please notice colors may vary slightly from those displayed here.

* Wheel rims are always black.



Invacare Australia Pty Ltd
 ABN 45 074 676 378
 PO Box 5002, 1 Lenton Place,
 North Rocks, NSW 2151,
 Australia.
 Freephone: 1800 460 460
 Fax: 02 8839 5353
 Email: orders@invacare.com.au
 Web: www.invacare.com.au

Invacare New Zealand
 PO Box 62-124, 4 Westfield Place,
 Mt. Wellington, Auckland,
 New Zealand.
 Freephone 0800 468 222
 Freefax 0800 807 788
 Email sales@invacare.co.nz
 Web: www.invacare.co.nz

Invacare Asia Ltd.
 PO Box 5002, 1 Lenton Place,
 North Rocks, NSW 2151,
 Australia.
 Phone: +61 2 8839 5330
 Fax: +61 2 8839 5331
 Email: asiasales@invacare.com
 Web: www.invacare.com



Yes, you can.

For more pre-sales information about this wheelchair, please consult the applicable website.
 For information about typical user description, intended use, adjustments and transportation, please refer to the user guide.